

# Lecture 15


## PD imaging and MRI magnets

Read section 5.14 from Smith and Webb for  
magnets

# Some facts about $T_1$ and $T_2$ weighing

## $T_1$ weighting:


- Fluids are dark (**why? - homework**), normal tissue is mid-grey or bright.
- Clear boundaries between the different tissues seen. **Anatomy scans.**



*$T_1$ -weighted*  
( $TR = 600$ ,  $TE = 11$ )

## $T_2$ weighting:

- Fluids are bright, normal tissue is dark.
- Abnormal collections of fluid are also bright. **Pathology scans.**



*$T_2$ -weighted*  
( $TR = 3800$ ,  $TE = 102$ )

# Proton density and contrast

# Proton density (PD)


- Proton density: Total number of excitable spins ( $N$ ) in the imaged volume.
- Decides the maximum signal that can be obtained from a particular tissue.  $N$  gives a measure of  $M_0$ .

Tissue	Proton density (a.u.)
Grey matter	85
White matter	70
CSF	100
Metastasis	85
Fat	100

# How will you choose $T_R$ and $T_E$ for proton density weighing?

- Choose very long  $TR$  (2500 ms) and very short  $TE$  (15 ms) to eliminate  $T_1$  and  $T_2$  dependency.
- Works comparatively better for fat and CSF (high PD. Why though?)

$$M_0 = \frac{N(\gamma\hbar)^2 B_0}{4kT}$$



# Compare how CSF and fat look in all three images



*T<sub>1</sub>-weighted*  
(TR = 600, TE = 11)

*Density-weighted*  
(TR = 3000, TE = 17)

*T<sub>2</sub>-weighted*  
(TR = 3800, TE = 102)

	T <sub>1</sub>	T <sub>2</sub>
CSF	2400	200
Fat	270	80

# Superconducting magnet ( $B_0$ )

- Must be spatially homogeneous. Why?
- Temporally stable (drift  $\sim 1\text{ppb}$ ) during MRI ( $\sim 40\text{ min}$ )



Typically **1.5T - 3T**


**Solenoid** (homogeneous field at the centre)


- For  $B_0 \sim 3\text{T}$ ,  $I = 350\text{A}$ . **Could the high current lead to a problem?**

Solution: Ti-Nb in a Cu matrix



# Cooling needed






- **Shielding:** reduce effect of 'stray field' outside magnet.
- Passive: put magnet in carbon-steel enclosure
- Active: shim coils

# RF coil ( $B_1$ )

- Same coil transmits and receives the signal.
- Highest efficiency by matching the size of the coil to the body part being imaged.



Knee coil




Head coil with mirror (for fMRI)

- **Body coil:** abdominal imaging

- We now understand why the steady (1.5T) and the RF magnets ( $\sim$  mT) are there.
- But why do we need the third magnet in the MRI set-up?

# What magnetic fields do we have?




- Steady magnetic field,  $B_0$  (initially aligns spins)
- RF magnetic field,  $B_1$  (excites aligned spins)
- Spatial information

# Gradient coils

Sections 5.8 and 5.9 from Smith and Webb

# Localizing the signal

We still don't have a way to distinguish the magnetic resonance signals coming from protons in different locations in the body. We need to localize the signal.



- Pulsed gradients  
(~ milliseconds)

- $G$  ( $\sim 10$  mT/m)  $\ll B_0$  ( $\sim 1.5$  T)
- Over 30 cm, total field variation is  $\sim 3$  mT (0.2% of  $B_0$ )
- Expressed in units of change in resonant frequency per cm (Hz/cm). For protons  $10$  mT/m = 4,258 Hz/cm.