# Lec 6: X-ray tube, attenuation photoelectric effect


# How does the x-ray spectrum change when the tube parameters are changed?

Think about what will happen to

- 1. Intensity
- 2. Average energy
- 3. Characteristic peaks


## Increase in tube voltage

#### Increase in tube voltage




#### Increase in tube current

#### Increase in tube current



### Increase in target (anode) material Z

#### Increase in target (anode) material Z




Relative Intensity

8


#### Increase in filter thickness

#### Increase in filter thickness



## Interaction of X-rays with tissues

#### Interaction with tissues leads to x-ray attenuation




Hendee, chapter 7
Smith and Webb, chapter 2

#### Attenuation mechanisms: absorption and scattering

Absorption: removal of x-ray photons from the beam


Scattering: change in direction of the photon, usually with reduced energy



#### Attenuation and contrast?

High contrast needs sufficiently high differential attenuation of x-rays in various tissues.

#### X-ray attenuation



How would you compare attenuation in two tissues?

$$\overline{I_{o}}$$

#### Linear and mass attenuation coefficients

Linear (
$$\mu$$
)  $\mu$ : 1/cm  
Mass ( $\mu_m = \mu/\rho$ )  $\mu_m$ : cm<sup>2</sup>/gm  


$$I(x) = I_o \exp(-\mu x)$$

$$I(x) = I_o \exp(-\mu x)$$

Watch out when solving problems!


#### Value of $\mu$ depends on x-ray energy

X-ray photon energy (keV)	Mass attenuation coefficient (cm <sup>2</sup> /g)		
	Compact bone	Muscle	
10	19.0	4.96	
15	5.89	1.36	
20	2.51	0.544	
30	0.743	0.154	
40	0.305	0.0677	
50	0.158	0.0409	
60	0.0979	0.0312	
80	0.0520	0.0255	
100	0.0386	0.0252	
150	0.0304	0.0276	
200	0.0302	0.0297	
300	0.0311	0.0317	




Which energy range should we choose for better contrast?

### Let's take a look at the y-axis scales...




# Plot x-ray intensity in bone and muscle as it progresses in x direction. Assume a tissue depth of 5cm.

Tissue	Density	Effective atomic	0.8-
	(g/cc)	number	0.6
Muscle	1.06	7.4	≤° 0.4
Fat	0.91	6.9	_ 0.4-
Bone	1.85	13.8	0.2-
	•		1




#### Recap: Value of μ depends on x-ray energy

X-ray photon energy (keV)	Mass attenuation coefficient (cm <sup>2</sup> /g)		
	Compact bone	Muscle	
10	19.0	4.96	
15	5.89	1.36	
20	2.51	0.544	
30	0.743	0.154	
40	0.305	0.0677	
50	0.158	0.0409	
60	0.0979	0.0312	
80	0.0520	0.0255	
100	0.0386	0.0252	
150	0.0304	0.0276	
200	0.0302	0.0297	
300	0.0311	0.0317	



Which energy range should we choose for better contrast?


#### Photoelectric effect




- Photon loses <u>all</u> its energy in one interaction with the tissue
- Watch out- the schematic diagram looks very similar to generation of characteristic x-rays.
- Inner shells

#### Auger electrons in photoelectric effect

# X-ray photon ejects electron



Transition of electron from higher to lower level emits another photon



Auger transition ejects another electron instead of photon

- Transition from higher (say, L) to lower energy shell (say, K) ejects another electron (usually from the same higher energy shell).
- Causes ionization of tissue.

K.E. 
$$_{Auger} = E_{Bi} - 2E_{Bo}$$