

Lec 10

(MRI: the role of B_0 and B_1 fields)

- A nucleus with either an odd number of protons or odd number of neutrons will have a net magnetic moment.
- Why does ^{14}N have a net magnetic moment then?

<i>Nuclide</i>	<i>Number of Protons</i>	<i>Number of Neutrons</i>
^1H	1	0
^2H	1	1
^{13}C	6	7
^{14}N	7	7
^{17}O	8	9
^{19}F	9	10
^{23}Na	11	12
^{31}P	15	16
^{39}K	19	20


Net bulk magnetization is along B_0

Bulk magnetization: $\mathbf{M} = \sum_{i=1}^N \mathbf{m}_i$


$$\langle M_z \rangle \neq 0, \langle M_x \rangle = 0, \langle M_y \rangle = 0$$

What happens during MRI? (1)


1. Hydrogen nuclei in tissue have “spin angular momentum” and associated magnetic moment.



2. In an external magnetic field (B_0), M_z lines up with B_0 (along z-axis).




3. A rotating magnetic field (B_1) pulse is applied along x-axis.




What happens during MRI? (2)

4. B_1 pulls away magnetization (M_z) from the z-axis with an angle α .




5. M_z rotates around z-axis at the “Larmor frequency”.




What happens in $B_1(t)$ field?

Magnetization in X-Y plane: RF field

- Apply RF pulse (10-100 MHz) B_1 along x-axis.
- B_1 ($\sim \mu\text{T}$ - mT) $\ll B_0$ ($\sim \text{T}$).




Net magnetization in xy plane (for $\pi/2$ pulse).




- M_z starts rotating about z-axis
- M_z returns to its equilibrium position along z-axis when RF field is turned off

If the magnitude of B_1 is much smaller ($\sim \text{mT}$) than B_0 ($\sim \text{T}$), then how is it possible that B_1 can flip some spins?

How do we detect M_z ?



Can't detect constant M_z
(with B_0 alone).





Can detect time-varying flux of
 M_z (generated using B_0 and B_1)

RF coils for generating B_1 can also detect the voltage
signal (**Faraday induction**).

$$V \sim - d\Phi/dt$$

We introduce a rotating reference frame. Why?





- In reality, each spin sees a slightly different magnetic field. Solving equations (that describe the time evolution of the magnetization vector, \mathbf{M}) becomes a nightmare.
- All MRI hardware has been designed to work assuming a rotating reference frame.

Frame rotates about z-axis with Larmor frequency.

Lab frame coordinates: x, y, z

Rotating frame coordinates: x', y', z'



What do we see if we are sitting on the rotating frame?

B_1 is constant.

B_0 vanishes.

