


# Lec 3: Generation of x-rays and X-ray tube

# Characteristic x-rays: have specific energies



1. Electron from cathode knocks out an inner shell (K) electron from the anode
2. Another electron from a higher energy shell (L) in anode fills the vacancy
3. Energy given up by  $L \rightarrow K$  transition is emitted as x-rays.


# Wavelength and energy of X-rays

Shell	Tungsten (keV)	Molybdenum (keV)
K	69.5	20
L	10.2-12.1	2.5-2.8
M	1.9-2.8	0.4-0.5


Example: Calculate the wavelength range of characteristic x-rays emitted during a transition from M level to K level in Tungsten.

If, E in keV and  $\lambda$  in nm, then  $E \text{ (keV)} = 1.24 / \lambda \text{ (nm)}$

# Nomenclature of characteristic lines




# X-ray production: Coolidge tube



Some of the kinetic energy of electrons hitting the target is converted into x-ray photons; the rest is dissipated as heat.

# Continuous energy x-rays: Bremsstrahlung



- Instead of knocking off an electron in the anode, electron from cathode changes its direction
- X-rays are emitted with continuous energies.


<http://physicsopenlab.org/2017/08/02/bremsstrahlung-radiation/>

$$E_{\text{x-ray}} = E_{\text{incident}} - E_{\text{final}}$$

# X-ray spectrum

# of photons  
in the beam

X-ray intensity



Diagnostic x-rays:  
20-150 keV